


NEW

FORESTADENT[®]
IMPLANTS 

OT EQUATOR

New abutment system
now available!



The cost-effective alternative.

One of the smallest attachment systems.

With a low vertical profile of 2.1 mm and diameter of 4.4 mm OT Equator is one of the smallest attachment systems on the market.











Female caps are retained by means of a stainless steel housing ranging in four levels of retention, making it easy to process at the dental laboratory or chairside in the dental office.

Your benefits at a glance.

- Cost-effective alternative to competitor products.
- Compatible for the Prowital implants.
- Compatible with Prowital screwdrivers.
- Easy to use.

OT Equator System for Implantats


Components:

Item	Description
	OT Equator Abutment Titanium + Tin Attachment; available in four gingival heights: 1, 2, 3, und 4 mm.
	Stainless Steel Housing
	Retentive Cap purple strong 2,7 kg
	Retentive Cap white standard 1,8 kg
	Retentive Cap pink soft 1,2 kg
	Retentive Cap yellow extra soft 0,6 kg
	Processing Cap
	Protective Disk
	Equator Impression Coping
	Equator Analog




Smart Box

The self-aligning retention cap

The new Smart Box is a retentive caps housing with innovative design. Thanks to a pivoting mechanism with a rotational core, it allows a passive insertion under extreme conditions, also up to 50° divergence.



Smart Box Components:

Item	Description
	Smartbox black Positioning Cap
	Pink Protective Disk
	Smart Box Housing with Black Positioning Cap


Available Kits

For OT Equator & Smart Box System

Order-No.	Description
1.278.35... ..10 = 1mm ¹ ..20 = 2mm ¹ ..30 = 3mm ¹ ..40 = 4mm ¹	OT Equator Abutment Kit - Ø 3,5 mm 1x OT Equator Titanium + Tin Attachment ¹ 1x Stainless Steel Housing 1x Protective Disk 4x Retentive Caps (1x strong, 1x standard, 1x soft, 1x extra soft)
1.278.43... ..10 = 1mm ¹ ..20 = 2mm ¹ ..30 = 3mm ¹ ..40 = 4mm ¹	OT Equator Abutment Kit - Ø 4,3 mm 1x OT Equator Titanium + Tin Attachment ¹ 1x Stainless Steel Housing 1x Protective Disk 4x Retentive Caps (1x strong, 1x standard, 1x soft, 1x extra soft)
1.278.0600	Equator Laboratory Set 1x Stainless Steel Housing 1x Processing Cap 1x Protective Disk 4x Retentive Caps (1x strong, 1x standard, 1x soft, 1x extra soft)
1.278.0605	Equator Laboratory Set Smart Box 1x Smart Box Housing with Black Positioning Cap 1x Pink Protective Disk 4x Retentive Caps (1x strong, 1x standard, 1x soft, 1x extra soft)

Auxiliaries

For retention caps

Order-No.	Item	Description
1.278.0010		Equator Instrument Multiuse Caps Inserter/Extractor Tool

1) Available in four gingival heights 1, 2, 3, and 4 mm.

FORESTADENT[®]
IMPLANTS 

DISTRIBUTOR:

FORESTADENT Implants GmbH
Westl. Karl-Friedrich-Straße 151
75172 Pforzheim / Germany

Phone: +49 (0)7231 / 459 - 250
service@forestadent-implants.com

RHEIN83

MANUFACTURER:

Rhein 83 srl • Via E. Zago, 10/ABC • 40128 Bologna (Italien)
Tel. (+39) 051 244 -510 / -396 • Fax (+39) 051 245238